



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Accessories For Bar Gel Polyacrylamide Electrophoresis Apparatus Circular

Product Image



Description

Accessories For Bar Gel Polyacrylamide Electrophoresis Apparatus Circular

Accessory 1

- Stabilized D.C. Power Supply : this is very ideal to accomplish for routine works coarse & fine control knobs are provided to set at any voltage from 0 to 300 & constant current up to 100 mA.

Accessory 2

- Stabilized D.C. Power Supply with digital control & display : This is provided with short circuit & over load protection with polarity reverse switch, coarse & fine

control knobs are provided to set at any voltage from 0 to 300 & constant current up to 100 mA.

Accessory 3

- Gel Slicer (Acrylic)

Accessory 4

- Tube Stand for Gel casting fixed

Accessory 5

- Levelling table with bubble eye

Accessory 6

- Glass Tube 10cm. x 5mm. I.D. set of 12

Accessory 7

- Glass Tube 10cm. x 6mm. I.D. set of 12

Accessory 8

- Silicon rubber grommets Set of 12

Accessory 9

- Rubber Bungs Set of 10

Accessory 10

- Polythene Clousers Set of 20

Accessory 11

- Syringe with Needle 1 Set

Accessory 12

- Gel destainer, diffusion type

Particulars	Catalog No.
Accessory 1	JA20780
Accessory 2	JA20781
Accessory 3	JA20782
Accessory 4	JA20783
Accessory 5	JA20784
Accessory 6	JA20785
Accessory 7	JA20786

Accessory 8	JA20787
Accessory 9	JA20788
Accessory 10	JA20789
Accessory 11	JA20790
Accessory 12	JA20791

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.