



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Accessories For Bar Gel Polyacrylamide Electrophoresis Apparatus Rectangular

Product Image



Description

Accessories For Bar Gel Polyacrylamide Electrophoresis Apparatus Rectangular

Accessory 1

- Stabilized D.C. Power Supply : this is very ideal to accomplish for routine works coarse & fine control knobs are provided to set at any voltage from 0 to 300 & constant current up to 100 mA.

Accessory 2

- Stabilized D.C. Power Supply with digital control & display : This is provided with short circuit & over load protection with polarity reverse switch, coarse & fine

control knobs are provided to set at any voltage from 0 to 300 & constant current up to 100 mA.

Accessory 3

- Silicon rubber grommets, Set of 12

Accessory 4

- Leveling table with bubble eye.

Accessory 5

- Gel cutter S.S. with acrylic handle

Accessory 6

- Tube stand for gel casting adjustable

Accessory 7

- Rubber Bungs Set of 10

Accessory 8

- Polythene Clousers Set of 20

Accessory 9

- Syringe with Needle 1 Set

Accessory 10

- Gel destainer, diffusion type

Particulars	Catalog No.
Accessory 1	JA20810
Accessory 2	JA20811
Accessory 3	JA20812
Accessory 4	JA20813
Accessory 5	JA20814
Accessory 6	JA20815
Accessory 7	JA20816
Accessory 8	JA20817
Accessory 9	JA20818
Accessory 10	JA20819

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.