



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Accessories For Submarine Gel Electrophoresis, Analyst Mechanism

Product Image



Description

Accessories For Submarine Gel Analyst Mechanism

Accessory 1

- Stabilized D.C. Power Supply : It is very ideal to accomplish for routine works, coarse & fine control knobs are provided to set at any voltage from 0 to 300 & constant current up to 100 mA.

Accessory 2

- Stabilized D.C. Power Supply with digital control & display : This is provided with short circuit & over load protection with polarity reverse switch. Coarse & fine

control knobs are provided to set at any voltage from 0 to 300 & constant current up to 100 mA.

Accessory 3

- Gel casting mechanism with leveling table.

Accessory 4

- Gel trays 13 x 13cm.

Accessory 5

- Gel trays 13 x 6 cm (2 nos.)

Accessory 6

- Gel trays 6 x 6.5cm.(4 nos.)

Accessory 7

- Acrylic Comb 15 well x 1.5mm

Accessory 8

- Acrylic Comb 8 well x 1.5mm

Accessory 9

- Acrylic Comb 2 wellx3mm

Accessory 10

- Platinum Electrode Assembly (2 nos.)

Accessory 11

- Developing Chamber (Black)

Particulars	Catalog No.
Accessory 1	JA20730
Accessory 2	JA20731
Accessory 3	JA20732
Accessory 4	JA20733
Accessory 5	JA20734
Accessory 6	JA20735
Accessory 7	JA20736
Accessory 8	JA20737
Accessory 9	JA20738

Accessory 10 JA20739

Accessory 11 JA20740

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.