



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Accessories For Submarine Gel Electrophoresis, Student Mechanism

Product Image



Description

Accessories For Submarine Gel Student Mechanism

Accessory 1

- Vari Pack : It is very simple to operate. It is very ideal for student practical classes & routine clinical work. it is designed to give 50, 100, 200 or 250 volts constantly indicated by LED on each point. Its maximum current is upto 100 mA.

Accessory 2

- Gel Casting mechanism with leveling table.

Accessory 3

- Gel trays 5x7cm (2 nos.)

Accessory 4

- Gel trays 10x7 cm.

Accessory 5

- Comb 8 well x 1.5mm

Accessory 6

- Comb 2 well x 3.0 mm

Accessory 7

- Comb stand

Accessory 8

- Developing Chamber for ethidium bromide staining

Accessory 9

- Gel Scoop U.V.T.

Particulars	Catalog No.
Accessory 1	JA20710
Accessory 2	JA20711
Accessory 3	JA20712
Accessory 4	JA20713
Accessory 5	JA20714
Accessory 6	JA20715

Accessory 7	JA20716
Accessory 8	JA20717
Accessory 9	JA20718

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.