



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

---

# Accessories For Submarine Gel Electrophoresis, Student Mechanism

## Product Image

---



## Description

---

### Accessories For Submarine Gel Student Mechanism

#### Accessory 1

- Vari Pack : It is very simple to operate. It is very ideal for student practical classes & routine clinical work. it is designed to give 50, 100, 200 or 250 volts constantly indicated by LED on each point. Its maximum current is upto 100 mA.

#### Accessory 2

- Gel Casting mechanism with leveling table.

#### Accessory 3

- Gel trays 5x7cm (2 nos.)

#### Accessory 4

- Gel trays 10x7 cm.

#### Accessory 5

- Comb 8 well x 1.5mm

#### Accessory 6

- Comb 2 well x 3.0 mm

#### Accessory 7

- Comb stand

#### Accessory 8

- Developing Chamber for ethidium bromide staining

#### Accessory 9

- Gel Scoop U.V.T.

<b>Particulars</b>	<b>Catalog No.</b>
Accessory 1	JA20710
Accessory 2	JA20711
Accessory 3	JA20712
Accessory 4	JA20713
Accessory 5	JA20714
Accessory 6	JA20715

Accessory 7	JA20716
Accessory 8	JA20717
Accessory 9	JA20718

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.