



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Anderson Pipette

Product Image



Description

Particle Size Determination Apparatus (Sedimentation pipette, Anderson) consisting of 500ml graduated cylinder, graduation from 0 to 20 cm and fitted with 10ml pipette with three way stopcock, the tip of pipette at the level of 0 mark at the cylinder.

-

Anderson Pipette

Catalog No. 100621

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.