



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Balance, Analytical, Mechanical Type

Product Image



Description

1. Digital (mechanical) indication with large figures for rapid and precise reading
2. Clearly-lit scale window the figures immediately above and below the measured value are also visible through green transparent shades.
3. Very short settling time and high weighing accuracy.
4. All controls conveniently located for efficient and reliable operation
5. Drop weight system with built-in drop weights in anti-magnetic stainless steel.
6. Zeroing device for setting the loaded machine to zero reading.
7. Stainless steel pan.
8. Weighing space with sliding side doors and fixed glass front.
9. Floor of weighing space also of glass.
10. Powder coat-painted pressure die cast aluminium housing.

Readability	0.01mg	0.1mg	0.1mg
Capacity	100mg	160mg	200mg
Optical Range	100mg	1000mg	1000mg
Scale Division	1mg	10mg	10mg
Built-in Mechanical Weight	0.1-99.9 gm	1-159 gm	1-199 gm
Accuracy of optical range	±0.01mg	±0.1 mg	±0.1 mg
Weighing principle	Substitution	Substitution	Substitution
Power Connection	220/230 volt AC	220/230 volt AC	220/230 volt AC

Catalog No.	Model
JW72501	Mechanical Analytical Balances, 100mg
JW72502	Mechanical Analytical Balances, 160mg
JW72503	Mechanical Analytical Balances, 200mg

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.