



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

---

# Ball and Ring apparatus, Brass

## Product Image

---



## Description

---

Ball And Ring Apparatus is used for demonstrating thermal expansion, comprising a captive brass ball of 19mm secured to a mounted brass ring by a chain. Ring mounted on rod with a wooden handle. The ball passes through the ring when cold but will not pass through after being heated.

### **Catalog No.    Particulars**

JZ00160            Ball and Ring apparatus, Brass

A quality product from JAPSON ideal for use in School Laboratories in developing countries of Asia, Africa and South America.

Standard specifications for school laboratory use. Ideal mix of Quality and Economy. Other configurations and superior quality available for specific use or for developed countries.

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.