



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Bar Gel Polyacrylamide Electrophoresis Apparatus Circular

Product Image



Description

Gel can be prepared in glass tubes. This mechanism can be used in controlled temperature for the analysis of heat liable proteins, isoenzymes, isoelectric focusing & other preparative studies. It is a traditional mechanism meant for running 8 gel tubes of 10cm. long. Acrylic pipes are used to make its circular shape cathode & anode electrode assembly is centrally mounted in circular shape, which provides uniform electrical flow. Three ventilators are provided in the lower buffer reservoir for cool air circulation. The elementary unit consists of anode reservoir, Cathode reservoir, Platinum Electrode assembly (Anode & Cathode)(2 Nos.), Glass Tubes 10 cm.(24 Nos.), Silicone Rubber Gromets Fixed & Spare(8+8 Nos.), Rubber Bungs (12 Nos.), Polythene Closers (24 Nos.), gel Casting stand, Syringe with Needle, Power Cord

Particulars

Catalog No.

Bar Gel Polyacrylamide Electrophoresis Apparatus
Circular

JA20770

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.