



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

# Barlows Wheel apparatus

## Product Image



## Description

---

Barlows Wheel Apparatus is used for demonstrating the conversion of electric energy into mechanical energy. Star-shaped copper disc rotatable vertically on pair of adjustable centers and dipping into mercury trough between the vertical limbs of a small horseshoe magnet on a wood base. Disc adjustable vertically on metal pillar which is connected to one of two terminals mounted on base.

The other terminal is connected to the mercury trough.

Supplied without Mercury.

Catalog No.	Particulars
JZ00200	Barlows Wheel apparatus

A quality product from JAPSON ideal for use in School Laboratories in developing countries of Asia, Africa and South America.

Standard specifications for school laboratory use. Ideal mix of Quality and Economy. Other configurations and superior quality available for specific use or for developed countries.

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.