



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

---

# Bell Experiment Apparatus (Bell in Bell Jar)

## Product Image

---



## Description

---

For demonstrating that sound cannot be transmitted in a vacuum. It consists of a 2 litre bell jar with a rubber stopper fitted with binding posts and connectors from which an electric bell is suspended. It operates on 3 to 6 volts, using either battery or a transformer. To be used with Air Pump (Vacuum Pump with Aluminium Plate)

Custom made configuration is available for developing countries or Ministry of Education Tenders

<b>Catalog No.</b>	<b>Particulars</b>
JP10365	Bell Experiment Apparatus (Bell in Bell Jar)

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.