



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Bernoulli Tube Set

Product Image



Description

A horizontal tube 500mm long with a bore of 10mm is constructed near its centre & short side tubes are provided in each section for the attachment of glass manometer tubes 270mm long with 8mm bore. The apparatus is supplied complete with three glass manometer tubes & rubber tubes for connecting manometer tubes.

Catalog No.	Particulars
-------------	-------------

JZ00270	Bernoulli Tube Set
---------	--------------------

A quality product from JAPSON ideal for use in School Laboratories in developing countries of Asia, Africa and South America.

Standard specifications for school laboratory use. Ideal mix of Quality and Economy. Other configurations and superior quality available for specific use or for developed countries.

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.