



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Bottle Top Dispenser, Premium Model

Product Image



Description

Premium Model :

All wetted parts are PTFE, FEP and Borosilicate Glass (3.3) only. This can be used for all laboratory reagents except HF.

General Product Information:

Description : Premium model

Piston : PTFE

Cylinder : Borosilicate Glass (3.3)

Volume adjustment Knob : Screw Types

Valve assembly : Glass ball & seat

Valve Housing : PTFE

Discharge assembly : PTFE

Delivery tube : FEP

Inlet tube : FEP

Calibration : Individually Calibrated and certified in lab easy calibration by user is also possible.

Accuracy & Reproducibility : In accordance with a ISO 8655 standards

Compatibility : Excellent Compatibility with all reagent except HF

Premium Model Available in 5 adjustable volume sizes :

(a) 0.25-2.5 ml **(b)** 0.5-5.0 ml **(c)** 1.0-10ml **(d)** 2.5-25.0 ml **(e)** 5.0-50.0ml

- Easy to use, simple construction and smooth, effortless plunger movement for bubble free dispensing.
- Piston : Specially designed PTFE Piston with Silicone O-ring ensures smooth, effortless piston movement and high accuracy.
- Each instrument comes with five adaptors to fit in most of the laboratory reagent bottles.

Standard threads : 30mm

Adaptors :28mm, 32mm, 36mm, 40mm, and 45mm.

- Nozzle cap prevents any unwanted drop on the work space.
- Easy to dis-assemble for cleaning and servicing.
- Glass Barrel is protected by a polypropylene sleeve. It also ensures no glass splitting in case of any accidental breakage.
- Valve Manifold : Specially designed, single piece PTFE valve manifold in the Premium Model ensures excellent chemical compatibility, leak proof valve functioning and excellent reproducibility.
- Instrument is fully autoclavable at 121 C, 15psi. No Need for disassembling the

instrument for autoclaving.

- Screw type volume adjustment knob with large pointer ensures precise, easy and quick volume adjustment.
- Calibration : Specially designed universal calibration tool is provided for convenient and quick In-lab user re-calibration. This is in compliance with GLO/ISO Norms.

Bottle Top Dispenser , Premium Model

Catalog No. 100224

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.