



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Bottles for Tissue Culture

Product Image



Description

Tissue Culture Bottles are ideal for large yield of cell growth and multiplication. These culture bottle are available either non-treated or surface treated which offers improved cell adhesion.

Catalog No.	Particulars
JW07871	2000ml, Surface Treated, Sterilized
JW07872	2000ml, Non-Treated, Sterilized

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.