




Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Capillary Tube Apparatus

Product Image



Description

Ideal for demonstrating the relationship between capillary pressure and the bore diameter of the capillary tube. The apparatus comprises a powder coated metal frame of height 90 mm holding six 150 mm capillary tubes of different bores vertically in a trough. The trough is filled with water and the difference in height of the resulting columns of water is immediately clear.

-

Capillary Tube Apparatus

Catalog No. 101155

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.