



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Charles Law Apparatus (Combined)

Product Image



Description

To demonstrate Charles Law i.e, the relationship between pressure and temperature of given mass of the gas at fixed volume. Also Known as Jollys Apparatus. Consist of Jollys air bulb connected to a glass mercury reservoir through flexible rubber tubing, mounting on a polished wooden stand with heavy cast metal base with leveling screw. To read the level directly a meter rule is fixed on a stand provided with a scale graduated 0-100cm x 1mm with the reading in the both directions. Mercury reservoir tube mounted on vertical support rod with the help of sliding bracket that can be positioned anywhere along the length of the rod. Supplied without mercury.

Custom made configuration is available for developing countries or Ministry of Education Tenders

Catalog No.	Particulars
JP10415	Charles Law Apparatus (Combined)

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.