



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Density Determination Kit for Solid and Liquid Samples

Product Image



Description

Density Kit for Solid and Liquid Test

Density Kit Suitable for all Weighing Balance, if weighing pan less than 125mm

If Weighing Balance not having Density function facility, Customer have to calculate manually.

According to the Archimedean Principle, Solid Body immersed in a liquid apparently loses as much of its own weight as the weight of the liquid it has displaced. This makes it possible to determine the unknown value. Depending on whether the body in quotient is liquid or solid.

The international system of units specifies Kg/m^3 as the unit of density. However, the unit g/cm^3 is better suited to lab purpose

STANDARD ACCESSORIES

- ◆ Beaker - 1No
- ◆ Stand for Beaker - 1No.
- ◆ Density Kit Bracket - 1No
- ◆ Hanging Holder - 1No.
- ◆ Crystal Ball - 1 No.
- ◆ Tweezer - 1 No.
- ◆ Thermometer - 1 No

Density Determination Kit for Solid and Liquid Samples
101137

Catalog No.

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.