



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Digital Melting Point Apparatus

Product Image



Description

Catalog No : 100123

Digital Melting Point Apparatus – Model 115

TECHNICAL SPECIFICATIONS

✓ DISPLAY	4 Digit Seven Segment LED display
✓ TEMPERATURE	0-300 °C
✓ RESOLUTION	1°C
✓ CAPILLARY GLASS TUBE	100mm long, 0.8 to 1.1mm internal diameter
✓ BOILING POINT	6x75mm Glass Tube
✓ SAMPLE SIZE	Fine powder properly dried at a temperature considerably below its melting point, 4 to 6 mm high
✓ OIL-BATH	Silicon oil bath beaker - 200 ml capacity with teflon cover fitted with oil immersion heater. It is having three holes for holding the capillaries, one for boiling point tube and one for thermometer
✓ STIRRER SPEED	Controllable with front knob
✓ POWER	230V \pm 10% AC, 50 Hz
✓ WEIGHT	6 Kg. (Approx.)

ACCESSORIES

- ✓ Box of Capillaries
- ✓ Glass Tubes for Liquid Sample (6 x 75 mm) - 2 Nos.
- ✓ Magnetic Capsule
- ✓ Bottle of Silicon Oil - 250 ml
- ✓ Glass Cylindrical Beaker with Teflon cover fitted with immersion heater
- ✓ Instruction Manual, Dust Cover

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.