



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Draining Tray

Product Image



Description

Material : Polypropylene

These autoclavable perforated trays can be placed in sink for washing, cleaning and draining of various laboratory glass and plastic apparatus.

Catalog No.	Particulars
101442	400x300x100mm

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.