



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

E/M by Helical Method Experiment Apparatus

Product Image



Description

Objective : To study charge of an electron.

Features : Kit Comprises of High Voltage Power Supply Unit with intensity, focus X, Y, deflection & Solenoid current controls.,br> Two metes are provided in Power Supply for Acceleration Voltage control & for solenoid current controls.

One Solenoid Unit Comprises of one Solenoid mounted on Wooden stand & CRT is kept in the solenoid pipe.

| Catalog No. | Particulars |
|--------------------|--|
| JE22110 | E/M by Helical Method Experiment Apparatus |

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.