



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

E/M by Thomson Method Experiment Apparatus

Product Image



Description

Objective : To study charge of an electron using thomson method.

Features : High voltage power supply with one meter provided to measure the deflection voltage. Complete with intensity, focus & deflection controls. It is provided with accessories such as compass box, set of Elnico magnets, teak wood stand with crt mounted, teak wood u stand with scale on both the arms for placing magnets.

Catalog No. Particulars

JE22100	E/M By Thomson Method Experiment Apparatus
---------	--

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.