



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Electro Convulsometer

Product Image



Description

Electro Convulsometer is used for applying maximal electro-shock through corneal electrodes provided. Instrument is used to study the anticonvulsant activity of phonation against maximal electro-sock induced convulsions in rat.

The instrument is provided with:

- Digital Voltmeter
- Analog Type Ammeter
- Analog Timer
- Multiplier
- Three pair of corneal Electrodes

-

Electro Convulsometer

Catalog No. 100644

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.