



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

E/m By Millikan's Oil Drop Method

Product Image



Description

The Principal elements of the assembly are specially designed condenser, a light source and a measuring microscope. All are mounted on a common support.

Consisting of a heavy cast iron base which is mounted on a steel pillar that carries the microscope and the lamp house. The microscope can be adjusted smoothly by slow motion screws in the horizontal as well as in the vertical direction. Atomizer is also fitted on the same carriage. The apparatus does not need any critical adjustment and the oil drops appear in microscope within seconds. Microscope is provided with micro meter scale.

Power Supply for Millikan's Oil Drop Apparatus:

Specially designed for use with Milliken's Oil Drop Apparatus. It give IC regulated D.C. output 0.300 volt continuously variable to be applied to condenser plates and a 6 volt at 3Amp. A.C. supply for the lamp house.

-

E/m By Millikan's Oil Drop Method

Catalog No. 101017

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.