



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Friability Test Apparatus, Digital, Single Drum

Product Image



Description

Programmable rotation time with digital display.
Motorized addition/removal of tablets without opening the drum.
Auto/Manual mode of operation.

Special low RPM motors less noise and constant speed even at very low or high voltage.
Auto stop with audio alarm.
Fully automatic.

Measures Obsession Strength of tablets with programmable timer. Rotation speed 25 rpm +/- 1, auto & manual mode of operation with 1 acrylic drums, Results are displayed on 3 digit LED display and it has 3 digit pre-stable control.

Model

Catalog No.

Friability Test Apparatus, Digital, Single Drum

JA20270

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.