



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

GM Counter Experimental Setup Experiment Apparatus

Product Image



Description

Instrument consists of power supply for g. M. Tube with digital readout meter. 6 digit counting system with its power supply. Controls for reset, start & stop through in panel. One presentable timer with 3 digit display & thumb wheel switch. Complete with g. M. Tube stand & one radio active source b.

Catalog No. Particulars

JE22120 GM Counter Experimental Setup Experiment Apparatus

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.