



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

GM Counter Experimental Setup Experiment Apparatus

Product Image



Description

Instrument consists of power supply for g. M. Tube with digital readout meter. 6 digit counting system with its power supply. Controls for reset, start & stop through in panel. One presentable timer with 3 digit display & thumb wheel switch. Complete with g. M. Tube stand & one radio active source b.

Catalog No. Particulars

JE22120	GM Counter Experimental Setup Experiment Apparatus
---------	--

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.