



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

---

# Hall Effect Experiment Apparatus

## Product Image

---



## Description

---

### Hall effect Experiment

**Features :** The set up consists of electromagnet of 7.5 Kg, Power supply for electromagnet 0-100V/4 Amp with meter. Germanium crystal mounted & its constant current power supply 0-20 mA fitted with two digital meters for Hall voltage & crystal current. Gauss meter digital 0-20 K. Gauss with Hall probe.

Catalog No.	Particulars
-------------	-------------

JE20200	Hall Effect Experiment Apparatus
---------	----------------------------------

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.