



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Hall Effect Experiment Apparatus

Product Image



Description

Hall effect Experiment

Features : The set up consists of electromagnet of 7.5 Kg, Power supply for electromagnet 0-100V/4 Amp with meter. Germanium crystal mounted & its constant current power supply 0-20 mA fitted with two digital meters for Hall voltage & crystal current. Gauss meter digital 0-20 K. Gauss with Hall probe.

Catalog No.	Particulars
JE20200	Hall Effect Experiment Apparatus

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.