



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Hares Apparatus on Stand

Product Image



Description

For comparing the densities of liquids by measuring the heights of two columns of liquid produced by applying suction to center limb. Comprising three limbed glass tube mounted on polished wooden stand with meter scale 45 cm long with 0 at the bottom. The scale is divided in centimeters and millimeters. Supplied without rubber tube and pinch cock. Complete with all accessories

Particulars

Catalog No.

Hares apparatus on stand

JZ01310

A quality product from JAPSON ideal for use in School Laboratories in developing countries of Asia, Africa and South America.

Standard specifications for school laboratory use. Ideal mix of Quality and Economy. Other configurations and superior quality available for specific use or for developed countries.

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.