



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Hot Plate Magnetic Stirrer, Euro

Product Image



Description

Salient Features

◆ MICROPROCESSOR CONTROL

The LMMS-300 hot plate magnetic stirrer is operated by a digital control CPU. The CPU regulates the heating up to 300 °C

◆ DIGITAL SPEED CONTROL

Stirring speed is also controlled by the CPU from the keypad while the stirring speed is indicated on the bar graph.

◆ PROGRAMABLE

Operator can setup the objective temperature, hold temperature, stirring time and speed.

◆ AUTOMATIC TEMPERATURE CONTROL

Connect the temperature probe and place into glass container being stirred to measure and maintain the object temperature.

Accuracy is ± 1 °C.

◆ STANDARD ACCESSORIES

Stir Bar, Temperature Probe, Instruction Manual

SPECIFICATIONS	MODEL	LMMS300
	Display	LCD
	Max. Volume (ml)	2000 ml
	Speed (RPM)	1250 rpm
	Top Plate Size	135 x 135 mm
	Plate Material	Steel
	Max.Temp. (Plate Surface)	300°C
	Temp. Accuracy	$\pm 1^{\circ}\text{C}(<100^{\circ}\text{C})\pm 3^{\circ}\text{C}(<300^{\circ}\text{C})$
	Timer	999 minutes
	Power supply	AC 220V / 50 Hz
	Dimensions LxWxH (mm)	230 x 180 x 120
	Weight	2.3 kg

Hot plate Magnetic Stirrer, Euro

Catalog No. 101122

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of

product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.