



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

# Hot Plate Magnetic Stirrer, Euro

## Product Image



## Description

### *Salient Features*

#### ◆ MICROPROCESSOR CONTROL

The LMMS-300 hot plate magnetic stirrer is operated by a digital control CPU. The CPU regulates the heating up to 300 °C

#### ◆ DIGITAL SPEED CONTROL

Stirring speed is also controlled by the CPU from the keypad while the stirring speed is indicated on the bar graph.

#### ◆ PROGRAMABLE

Operator can setup the objective temperature, hold temperature, stirring time and speed.

#### ◆ AUTOMATIC TEMPERATURE CONTROL

Connect the temperature probe and place into glass container being stirred to measure and maintain the object temperature.

Accuracy is  $\pm 1$  °C.

#### ◆ STANDARD ACCESSORIES

Stir Bar, Temperature Probe, Instruction Manual

MODEL	LMMS300
Display	LCD
Max. Volume (ml)	2000 ml
Speed (RPM)	1250 rpm
Top Plate Size	135 x 135 mm
Plate Material	Steel
Max.Temp. (Plate Surface)	300°C
Temp. Accuracy	$\pm 1^{\circ}\text{C}(<100^{\circ}\text{C})\pm 3^{\circ}\text{C}(<300^{\circ}\text{C})$
Timer	999 minutes
Power supply	AC 220V / 50 Hz
Dimensions LxWxH (mm)	230 x 180 x 120
Weight	2.3 kg

Hot plate Magnetic Stirrer, Euro

Catalog No. 101122

## Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of

product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.