



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Inverted Tissue Culture Microscope, Superior

Product Image



Description

Body:

The binocular observation tube assembled from coated prisms is 30 degree inclined, with eyepiece position 280mm raised above the desk surface for greater comfort to the observer with qadruple revolving turret nosepiece with positive click stops.

Stage:

1. Oversize plain stage 180x150mm with two extension plates 150mm x 65mm for large bottles & microsoft work. Two stage clips are provided in drilled holes to accommodate specimen upto 40mm. A removable central insert gives 65 mm dia stage aperture.
2. Co-axial, cross movement mechanical stage 180x150mm.
3. Low positioned co-axial mechanical stage 180x50mm with three different aperture plates.

Focussing

Coarse Focussing:

Two large focus knobs located low down at a convenient position assure accurate focusing. The back-lash free rack & pinion mechanism has 18mm range.

Fine Focussing:

Highly sensitive based on screw-lever mechanism reading to 0.002mm.

Illumination:

Based on Koehler system is provided by a post-mounted high intensity high intensity 6 volts 20 watts Halogen illuminator, powered from base built in variable intensity controlled transformer. The lamp house incorporating a focusable condenser lens with iris diaphragm & filter-carrier is mounted on a pole on which it can be positioned up or down to accommodate any type of culture vessel in a vertical range of 140 mm working distance lens is 20mm

Eyepiece:

HWF 10X (paired) with eye-guards, for genuine wide-field observation.

Objective:

Achromatic 5X, 10X, LWD 20X & LWD 40X. Phase contrast specially designed perfect flat-field objectives with long working distance corrected for bottom thickness of bottles, flasks or chambers

Box:

Supplied in thermocol box.

Model

Catalog No.

Inverted Tissue Culture Microscope,
Superior

JM005900

Note: Available with Trinocular head

Note: Choice of two mechanical stages are also available as optional.

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.