



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

# Jaeger's Surface Tension Apparatus

## Product Image



## Description

---

Having a wide mouth reagent bottle fitted with a cork having two holes, in one hole a dropping funnel is fitted. Complete with manometer, capillary tube. All mounted on suitable wooden stand. Without aspirator.

Jaeger's Surface Tension Apparatus

**Catalog No.** 100410

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.