



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Laboratory Stirrer, High Torque, Geared Motor

Product Image



Description

Safe and a versatile unit. A very high torque is developed at a relatively low operating speed. It is suitable for affectively stirring high viscosity fluids and pastes. Geared Motor 1/20 H.P., 4000 rpm, continuous type motor is coupled to speed reducing gear unit giving output speed of 150 rpm at 220 volts AC/DC. Stirrer Shaft 8 mm. dia., 350 mm long shaft fitted with 4 bladed turbine impeller made of stainless steel. Speed Control Variable speed controller, mounted in a metal housing. Digital Speed controller also available. Stirring Capacity 2 to 3 kgs. of thick materials.

Catalog No.Particulars

JA23140	Laboratory Stirrer, High Torque, Geared Motor
JA23141	Laboratory Stirrer, High Torque, Geared Motor, Digital RPM Indicator & Controller From 200 To 2500 RPM

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.