

Address: JAMBU PERSHAD & SONS

6275/22 Nicholson Road, Ambala Cantt, Haryana, INDIA

Pin: 133001

Email:

sales@japson.com japsonambala@yahoo.com **Website:** www.japson.com

Phone:

+91-171-4006897

# **Laboratory Tray, Plastic**

## **Product Image**



## **Description**

#### Material: Polypropylene

Needless to say tray moulded in polypropylene is a multipurpose tray that can be used for sterilizing drying glassware porcelainware and other laboratory instruments. These steam autoclavable trays have tapered walls which prove good nesting.

Catalog No.	Size (mm)
JW11801	450x350x75
JW11802	375x300x75

### **Disclaimer**

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.