



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

---

# Laser Experimental Set Up (He-Ne Laser) Experiment Apparatus

## Product Image

---



## Description

---

Objective : Study the wavelength & angular divergence of beam.

Set-up consists of following apparatus :

He-Ne Laser Source with Power Supply

Set of 13 slides

Graticule

Grating 500 lines.

Optical Bench 1.5 meter.

Graph Paper with screen & scale.

Convex Lens 2 pcs. (FL 10 cm & 20 cm).

**Catalog No.    Particulars**

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.