



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Linear Expansion Apparatus

Product Image



Description

This Linear Expansion testing device makes studying the expansion of metals simple and instantly visible to the student. A metal rod is inserted into the steam chamber and steam is applied. Expansion of the rod is immediately registered by the indicator that is mounted on one end. The heat jacket has three inlet tubes for water intake, outlet and thermometer insertion. Steam source is not provided.

Particulars	Catalog No.
Linear Expansion Apparatus	JZ01780

A quality product from JAPSON ideal for use in School Laboratories in developing countries of Asia, Africa and South America.

Standard specifications for school laboratory use. Ideal mix of Quality and Economy. Other configurations and superior quality available for specific use or for developed countries.

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.