



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Magnetic Needle

Product Image



Description

PTFE/Teflon Coated rotor/needle

Strong chemical-resistance

Acid and alkali resistant

Insulated and low co-efficient

To be used with Magnetic Stirrer and Magnetic Stirrer with Hot Plate.

Particulars	Catalog No.
Magnetic Needle on Stand	JZ01930
Magnetic Needle only	JZ01931

A quality product from JAPSON ideal for use in School Laboratories in developing countries of Asia, Africa and South America.

Standard specifications for school laboratory use. Ideal mix of Quality and Economy.

Other configurations and superior quality available for specific use or for developed countries.

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.