



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Mammalian Heart Perfusion Assembly

Product Image



Description

A versatile assembly serving both for heat perfusion and as an all purpose thermostatic Organ Bath. The transparent Perspex bath is 12"x9"x8" and rests on a metal stand: a 1"x3/8" thick edge at the top takes clamps etc carrying the uprights and accessories. Temperature of the bath is regulated (CONSTANT). Thermostatic sensitivity $\pm \frac{1}{2}^{\circ}$ C. The stirrer unit is completely independent, and dips into the bath in one corner, and fixes on its special stainless steel upright. A Perspex connection box fitted to the bath frame, and carries all the electric leads and connections. Two pilot lamps indicate the thermostat and stirrer operation.

Other accessories for heart perfusion assembly.

Extra accessories to enable the unit to be used as all purpose Organ Bath (2 unit) are:-

- (a) Double set of glass inner vessels.
- (b) Oxygen tubes.

- (c) Warming coils.
- (d) Two stainless steel capillary levers.
- (e) Four stainless steel uprights with boss heads.
- (f) X Block

-

Mammalian Heart Prefusion Assembly

Catalog No. 100658

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.