



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

---

# Mayer's Viscosity Apparatus

## Product Image

---



## Description

---

In this, we can determine viscosity of liquid by oscillating disc method described by MAYER. A flat disc is suspended horizontally by a phosphorous bronze suspension which is attached to a rod rigidly fastened to the center of the disc. This rod carries a cross bar whose ends have a screw thread, along which two masses are screwed to balance the disc horizontally. A small concave mirror is fastened to the rigid support, complete mounted on base, but without lamp and scale outfit.

Mayer's Viscosity Apparatus

**Catalog No. 100412**

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change

them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.