



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

---

# Melting Point Apparatus, Digital, Fully Automatic

## Product Image

---



## Description

---

packed & easy to use, this is an indispensable bench top tool. This equipment has an arrangement to hold the capillary containing the sample in uniformly heating media. The temperature is indicated on a digital display readable in 1.0°C. At the point, when the sample melts, the heating is automatically cut-off and the digital display showing the temperature remains on the indicator till it is cleared. No Silicon oil or other bath is required as heating media. Range 2.0°C above room temperature to 350°C.

Fully automatic.  
Latest technology.  
Range from 40°C to 350oC.  
Highly accurate.  
Fully automatic.

### Technical Specification:

<b>Display:</b>	Digital 3½ digit, Bright LED
<b>Accuracy:</b>	1%
<b>Operation:</b>	Fully automatic
<b>Dimensions:</b>	405x360x235mm (approx)
<b>Weight:</b>	5Kg (approx.)
<b>Power:</b>	230V ± 10%AC, 50Hz

### Accessories:

Capillary Tube - 100 No.  
Cleaning Tool - 1 No.

### Model

Melting Point Apparatus, Digital, Fully Automatic

### Catalog No.

JA90150

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.