



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Melting Point Apparatus With Cooling System

Product Image



Description

A single knob on the front panel incorporate the mains On/off switch this gives a temperature range from ambient to 350°C with an arbitrary scale of 0 to 10. A rapid heat switch allows higher temperature to be obtained in a shorter time. A lamp provides uniform shadow less illumination of the sample tubes, which are viewed through a 25 mm diameter-

magnifying lens. A white or dark background can be selected with a control at the rear of the apparatus.

Construction:

An aluminum top casting enclosed the heating block & minimize the transfer of heat. The rest of the case is made from steel, covered with an epoxy/polyester paint finish, which is highly durable for rapid cooling of the apparatus between tests. A water cooled plug passes into a hole in the heating block, which provides rapid cooling for the next sample.

Particulars

Catalog No.

Melting Point Apparatus With Cooling System

JA90165

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.