



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Overflow Vessel (Eureka Can), Metal

Product Image



Description

G.I Sheet hammer-tone finish,the top edge is flared to add strength.
Used in Archimedes experiment.
Available in different metals on demand.

Custom made configuration is available for developing countries or Ministry of Education
Tenders

Catalog No.	Particulars
JP10310	Overflow Vessel, Metal

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.