

Address: JAMBU PERSHAD & SONS

6275/22 Nicholson Road, Ambala Cantt, Haryana, INDIA

Pin: 133001

Email:

sales@japson.com japsonambala@yahoo.com **Website:** www.japson.com

Phone:

+91-171-4006897

PCR Thermal Cycler, Digital

Product Image



Description

Catalog No : 100134

PCR Thermal Cycler, Digital - Model 240

TECHNICAL SPECIFICATIONS

✓ Block capacity

✓ Temp range

✓ Max Heating rate

✓ Max cooling rate

✓ Temperature display accuracy

✓ Accuracy

✓ Accuracy

64 x 0.2m1

4°C — 99°C

< 2.5°C/s

< 2.5°C/s

< 2.5°C/s

< + 2.5°C/s

✓ Accuracy
✓ Temperature fluctuation
< ± 2.5°C/s</p>
< ± 0.1 °C</p>

✓ Uniformity 95 °C h < \pm 0.5 °C (20 seconds later)

 $20 \,^{\circ}\text{C} - 72 \,^{\circ}\text{C} \, \text{h} < \pm 0.3 \,^{\circ}\text{C}$

(20 Seconds Later)

✓ Max program unit
 ✓ Max Cycle steps
 ✓ Max Cycle No
 ✓ Max Constant temperature
 ✓ Stored program No
 ✓ Stored program No
 ✓ Stored program No

✓ Dimensions 370mm (L) x 249mm (W) x I 80mm (H)

✓ Weight 5 kg

ACCESSORIES

- ✓ Adjustable pressure hot lid,
- ✓ Heat lid with pressure alarm device,
- Convenient and flexible module
- ✓ Innovative module wire socket design achieves module replacement without wire.
- ✓ The unique left-right design for amplification area and operating area.
- Coloured and Wide Screen displaying Multiple details.

Perform 20000 times non-breakdown operation.

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.