



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Perimeter, Lister type

Product Image



Description

The instrument automatically records field to 90 ° in all meridians. The crier mechanism is smooth, silent and invisible. The patient perceives only the test object to move. Spring check, giving an audible click to each 15° of rotation of the arc in either direction which renders close scrutiny of dial plate unnecessary, and ensures accurate registration. The device can be thrown out of action if silent movement is desired. The marker a needle is

housed in a retractable tube, sprung to sheath the pint when actually marking. A set of interchangeable and permanent washable test objects is included with the instrument. The objects subtend angles of 4° , 1° , $\frac{1}{2}^\circ$, $\frac{1}{4}^\circ$. Supplies complete with transformer to work on 220 volt A/C mains with 100 charts.

-

Perimeter, Lister type

Catalog No. 100660

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.