




Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Pipette Tray

Product Image



Description

Material : PS

These pipette trays can be used in drawers or on the bench top for storing pipettes and can also be used for soaking them before cleaning. Can accommodates 7 pipettes with overall stem diameter not exceeding 20 mm longitudinally or 16 pipettes with stem diameter not exceeding 10 mm laterally while the center support holds 12 short pipettes in a perpendicular direction. Pipettes can easily be places in and taken out from the troughs in the tray even with gloves.

-

Catalog No.	Particulars
101471	Pipette Tray, 283X216X40mm

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.