



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Platform Gel Rocker

Product Image



Description

Gel rocker can be used for many applications like gel staining, culture aeration, blotting membrane etc. A rubber mat is provided for no slipping surface for the platform. The movement is controlled very precisely with the D.C. Motor, with 0 - 60 minute timer for auto off within a time limit required. Gel rocker is provided with 30 x 30cm Platform size and having see saw type shaking system.

Particulars**Catalog No.**

Platform Gel Rocker

JA21000

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.