




Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001


Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Polensky Apparatus

Product Image





Description

Material:-Borosilicate Glass

**Catalog
No.**

Particulars

JG28401	Rubber cork type, for the determination of butter fat and analysis of oil, fats, waxes with rubber cork fitting with clamp, stand and burner.
JG28402	Joint type, for the determination of butter fat and analysis of oil, fats, waxes with standard joint fitting with clamp, stand and burner.

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.