



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Ria Vial

Product Image



Description

Material: Polypropylene/Polystyrene

These Polypropylene / Polystyrene tubes are generally used in RIA Coagulation & bacteriology. The Polypropylene tubes are clear and can withstand 3000 RCF during centrifugation while the transparent Polystyrene tubes can withstand 1400 RCF. The PP tubes can accept common acids, alkalies & solvents at room temperature while the PS tubes can accept mild bases, weak acids but not organic solvents.

Catalog No.	Particulars
JW06101	12x75mm in PS
JW06102	12x75mm in PP
JW06103	Cap. For 12mm RIA Tube

Disclaimer

The Product details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.