



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Rotary Microtome, Erma Type

Product Image



Description

It is a standard Erma Type Microtome. It has a well designed mechanism with latest cam drive system, enclosed in a metal cover and fitted on a heavy base.

1. Feed adjustment : 1 to 25 microns in step of 1 micron
2. Feed excursion : 25 mm.

3. Opening of object clamp : 30x15mm.
4. Overall size : 35x27x22 cms. approx.

Particulars**Catalog No.**

Rotary Microtome, Erma Type

JM00440

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.